



## Aesthetics - MUSTER

### LIGHT4YOU

Quick spot 60000 PULSED LIGHT-painless-a breeze



**muster**

✉ bbazar@globelife.com

☎ +39 0331 1706328



### LIGHT4YOU

Quick spot 60000 PULSED LIGHT-painless-a breeze

- 7 touch screen color display "
- 60000 spot
- Large size of the treated area for single Spotlight: ben 7.5 cm<sup>2</sup>
- Combined cooling air + water and thermoelectric cooler
- Treatment: the cooling rate of the handpiece and its sophisticated technology, guarantee a significant reduction of the treatment times.

#### PROFESSIONAL AIR And WATER COOLED EQUIPMENT

Pulsed light, an innovative, fast and effective, non-invasive hair removal treatment. Using the energy focused Xenon light, the application acts directly on the hair bulb with amazing results. The high intensity photo-epilation IPL exploits the principle of selective photothermolysis, generating a series of flashes that inhibit the cells responsible for hair development, contributing to the progressive and permanent removal of unwanted hair.

#### COMPUTERIZED MEDICAL HISTORY

Computerized medical history is simple and intuitive. Just tap the video at the desired function and the device allows you to easily consider all parameters to set proper treatment — all in a matter of seconds. LIGHT4YOU determines the appropriate level of power and the specific filter to use for treatment, by selecting, in a few steps a video, its parameters such as skin color and hair color and appearance. At the end of the video will display, for your appropriate verification, a summary of selected data and will proceed to treatment.

#### PHOTO-REJUVENATION TREATMENT

(optional)-not to be used by the operator in Italy (D.M. 110/2011).

The photo-rejuvenation is achieved thanks to the beneficial effects produced by pulsed light treatment, which stimulate the affected areas, awakening their proper and natural functioning. The innovative technologies offered by LIGHT4YOU showed visible improvements in skin abnormalities such as: rosacea, wrinkles, superficial capillaries, hypotonicity of the epidermis. The bright flash is able to activate the skin, stimulating the

production of new collagen, which revitalizes the skin texture and consistency.

- Spot 15 mm x 50 mm (7.5 cm<sup>2</sup>)
- Filters (2)
- Photo-epilation:
  - 600 nm filter (suitable for skin types I-II-III)
  - 640 nm filter (fits phototype IV-V-VI)
- Energy emitted: adjustable from 0.5 to 18 J/cm<sup>2</sup>
- Amplitude spot 20 ms
- Single or multiple spot spot Function (continuous output commercials to speed up the processing)
- Pause between the commercials in multiple spot mode: adjustable from 1 s to 4 s
- Charging time after the pulse: 250 ms (0.25 seconds)
- Light unit cooling system: combined water + air + thermoelectric cooler
- Crystal temperature setting of emission skin contact: -5 ° C to 5° C
- Lifetime 60000 spot lamp
- Average power consumption in operation: 180 W
- Maximum power charging: 2200 W
- Voltage 230 V/50 Hz